




TIFA(Phospho-Thr9) rabbit pAb

Catalog No	YP-Ab-10483
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TIFA T2BP
Immunogen	Synthesized peptide derived from human TIFA(Phospho-Thr9)
Specificity	This antibody detects endogenous levels of TIFA(Phospho-Thr9) at Human, Mouse, Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.138% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T2BP; T6BP; TIFA; TIFAA; TRAF2 binding protein; TIFA T2BP; TIFA(Phospho-Thr9); TRAF-interacting protein with FHA domain-containing protein A (Putative MAPK-activating protein PM14) (Putative NF-kappa-B-activating protein 20) (TRAF2-binding protein)
Observed Band	20kD
Calculated Molecular Weight	20kD
Cell Pathway	Cytoplasm . Colocalizes with lysosomal marker LAMP2 following homooligomerization and subsequent activation. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of cell lysis treated or untreated by LPS 100ng/mL 30min, using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000