



IRAK4(p-T345/S346) rabbit pAb

Catalog No	YP-Ab-10484
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	IRAK4
Immunogen	Synthesized peptide derived from human IRAK4(p-T345/S346)
Specificity	This antibody detects endogenous levels of IRAK4(p-T345/S346) at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.140% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IPD1;IRAK 4;IRAK4;NY REN 64;REN64;IRAK4(p-T345/S346);Interleukin-1 receptor-associated kinase 4 (IRAK-4) (EC 2.7.11.1) (Renal carcinoma antigen NY-REN-64)
Observed Band	52 kDa
Calculated Molecular Weight	52 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
