



QIK/QSK (Phospho-Thr175/163) rabbit pAb

Catalog No	YP-Ab-10487
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	SIK1 SIK SNF1LK
Immunogen	Synthesized peptide derived from human QIK/QSK (Phospho-Thr175/163)
Specificity	This antibody detects endogenous levels of QIK/QSK (Phospho-Thr175/163) at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.143% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SIK 2;salt inducible kinase 2;Qin induced kinase;QIK;LOH11CR11;QIK/QSK; SIK1 SIK SNF1LK;QIK/QSK (Phospho-Thr175/163);Serine/threonine-protein kinase SIK1 (EC 2.7.11.1) (Salt-inducible kinase 1) (SIK-1) (Serine/threonine-protein kinase SNF1-like kinase 1) (Serine/threonine-protein kinase SNF1LK)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Locates to cytoplasm when inactive following cAMP-induced phosphorylation, probably by PKA (PubMed:29211348).
Tissue Specificity	Fetal lung,Spleen,Testis,
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein., cofactor:Magnesium.,enzyme regulation:Activated by phosphorylation on Thr-182 by STK11 in complex with STE20-related adapter-alpha (STRAD alpha) pseudo kinase and CAB39.,function:Transient role during the earliest stages of myocardial cell differentiation and/or primitive chamber formation and may also be important for the earliest stages of skeletal muscle growth and/or differentiation. Potential role in G2/M cell cycle regulation. Inhibits CREB activity by phosphorylating and repressing the CREB-specific coactivators, CRT1-3.,



PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein
kinase family. AMPK subfamily.,similarity:Contains 1 protein kinase domain.,
similarity:Contains 1 UBA domain.,subcellular location:Translocates to the
cytoplasm on phosph

Background


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of 1 HeLa cell, 2 LPS
100ng/mL 30min treated ,using primary
antibody at 1:1000 dilution. Secondary
antibody(catalog#:RS23920) was diluted at 1:10000