



USP28 (Phospho-Ser714) rabbit pAb

Catalog No	YP-Ab-10526
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	USP28 KIAA1515
Immunogen	Synthesized peptide derived from human USP28 (Phospho-Ser714)
Specificity	This antibody detects endogenous levels of USP28 (Phospho-Ser714) at Human, Mouse, Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.182% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Deubiquitinating enzyme 28; EC:3.4.19.12; KIAA1515; Ubiquitin carboxyl-terminal hydrolase 28; Ubiquitin thiolesterase 28; USP28; USP28 KIAA1515; USP28 (Phospho-Ser714); Ubiquitin carboxyl-terminal hydrolase 28 (EC 3.4.19.12) (Deubiquitinating enzyme 28) (Ubiquitin thiolesterase 28) (Ubiquitin-specific-processing protease 28)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleoplasm .
Tissue Specificity	Brain, Eye,
Function	catalytic activity: Ubiquitin C-terminal thioester + H ₂ O = ubiquitin + a thiol., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Belongs to the peptidase C19 family., similarity: Contains 1 UIM (ubiquitin-interacting motif) repeat.,
Background	The protein encoded by this gene is a deubiquitinase involved in the DNA damage pathway and DNA damage-induced apoptosis. Overexpression of this gene is seen in several cancers. [provided by RefSeq, Oct 2016],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various,using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000