



MST3 (Phospho-Thr18) rabbit pAb

Catalog No	YP-Ab-10547
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	STK24 MST3 STK3
Immunogen	Synthesized peptide derived from human MST3 (Phospho-Thr18)
Specificity	This antibody detects endogenous levels of MST3 (Phospho-Thr18) at Human, Mouse, Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.203% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MST3; STK24 MST3 STK3; MST3 (Phospho-Thr18); Serine/threonine-protein kinase 24 (EC 2.7.11.1) (Mammalian STE20-like protein kinase 3) (MST-3) (STE20-like kinase MST3) [Cleaved into: Serine/threonine-protein kinase 24 36 kDa subunit (Mammalian STE20-like protein kinase 3 N-terminal) (MST3/N); Serine/threonine-protein kinase 24 12 kDa subunit (Mammalian STE20-like protein kinase 3 C-terminal) (MST3/C)]
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	Cytoplasm. Nucleus. Membrane. The truncated form (MST3/N) translocates to the nucleus. Colocalizes with STK38L in the membrane.
Tissue Specificity	Isoform A is ubiquitous. Isoform B is expressed in brain with high expression in hippocampus and cerebral cortex.
Function	catalytic activity: ATP + a protein = ADP + a phosphoprotein., cofactor: Manganese. Only isoform A requires manganese., function: Protein kinase that act on both serine and threonine residues., PTM: Autophosphorylated. Isoform B is activated by phosphorylation by PKA., similarity: Belongs to the protein kinase superfamily. STE Ser/Thr protein kinase family. STE20 subfamily., similarity: Contains 1 protein kinase domain., tissue specificity: Isoform A is



ubiquitous. Isoform B is expressed in brain with high expression in hippocampus and cerebral cortex.,

Background

This gene encodes a serine/threonine protein kinase that functions upstream of mitogen-activated protein kinase (MAPK) signaling. The encoded protein is cleaved into two chains by caspases; the N-terminal fragment (MST3/N) translocates to the nucleus and promotes programmed cells death. There is a pseudogene for this gene on chromosome X. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2013],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various,using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000