



Tryptase-3 Polyclonal Antibody

Catalog No	YP-Ab-10597
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	TPSD1
Immunogen	Synthesized peptide derived from Tryptase-3 at AA range: 51-100
Specificity	Tryptase-3 Polyclonal Antibody detects endogenous levels of Tryptase-3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000, ELISA 1:10000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Delta tryptase;HmMCP 3 like tryptase III;Mast cell mMCP 7 like;MCP7 LIKE; MCP7L1;MMCP 7L;TPSD1;Tryptase 3;Tryptase delta;tryptase delta 1;Tryptase-3;Tryptase delta (EC 3.4.21.59) (Delta-tryptase) (HmMCP-3-like tryptase III) (Mast cell mMCP-7-like) (Tryptase-3)
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	Secreted . Released from the secretory granules upon mast cell activation. .
Tissue Specificity	Expressed in colon, lung, heart and synovial tissue. May be specific to mast cells.
Function	catalytic activity:Preferential cleavage: Arg-I-Xaa, Lys-I-Xaa, but with more restricted specificity than trypsin.;caution:Although PubMed:11174199 reported this as a pseudogene, PubMed:12391231 showed it is expressed and has proteolytic activity when expressed in bacterial cells.;function:Tryptase is the major neutral protease present in mast cells and is secreted upon the coupled activation-degranulation response of this cell type.;similarity:Belongs to the peptidase S1 family. Tryptase subfamily.;similarity:Contains 1 peptidase S1 domain.;subcellular location:Released from the secretory granules upon mast cell activation.;subunit:Homotetramer.;tissue specificity:Expressed in colon, lung,



heart and synovial tissue. May be specific to mast cells.,

Background

tryptase delta 1(TPSD1) Homo sapiens Tryptases comprise a family of trypsin-like serine proteases, the peptidase family S1. Tryptases are enzymatically active only as heparin-stabilized tetramers, and they are resistant to all known endogenous proteinase inhibitors. Several tryptase genes are clustered on chromosome 16p13.3. These genes are characterized by several distinct features. They have a highly conserved 3' UTR and contain tandem repeat sequences at the 5' flank and 3' UTR which are thought to play a role in regulation of the mRNA stability. Although this gene may be an exception, most of the tryptase genes have an intron immediately upstream of the initiator Met codon, which separates the site of transcription initiation from protein coding sequence. This feature is characteristic of tryptases but is unusual in other genes. Tryptases have been implicated as mediators in the pathogenesis of asthma and other allergic and

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images




Western Blot analysis of HEPG2 cells using Tryptase-3 Polyclonal Antibody diluted at 1:800. Secondary antibody(catalog#:RS0002) was diluted at 1:20000




Immunohistochemical analysis of paraffin-embedded human-tonsil, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-tonsil, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-skin, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-skin, antibody was diluted at 1:200