




Pan glycine myristoylation Polyclonal Antibody

Catalog No	YP-Ab-10633
Isotype	IgG
Reactivity	Species independent
Applications	WB;ELISA
Gene Name	
Immunogen	Synthetic myristoylated peptide
Specificity	The antibody detects endogenous glycine myristoylated protein
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000, ELISA 1:10000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Pan glycine myristoylation
Observed Band	130 80 35kD
Calculated Molecular Weight	130 80 35kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of KB, MOUSE-LUNG cells using Antibody diluted at 1000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000