



Tetranectin Polyclonal Antibody

Catalog No	YP-Ab-10709
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC;IF;ELISA
Gene Name	CLEC3B TNA
Immunogen	Synthetic peptide from human protein at AA range: 11-60
Specificity	The antibody detects endogenous Tetranectin
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC-p 1:50-200, ELISA 1:10000-20000. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CLEC3B;TNA;TN;Plasminogen kringle 4-binding protein;C-type lectin domain family 3 member B;Tetranectin;CLEC3B TNA;Tetranectin (TN) (C-type lectin domain family 3 member B) (Plasminogen kringle 4-binding protein);Tetranectin (TN;Plasminogen kringle 4-binding protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted.
Tissue Specificity	Found in plasma.
Function	function:Tetranectin binds to plasminogen and to isolated kringle 4. May be involved in the packaging of molecules destined for exocytosis.,mass spectrometry: PubMed:10614823,online information:Tetranectin, similarity:Contains 1 C-type lectin domain.,subunit:Homotrimer.,tissue specificity:Found in plasma.,
Background	function:Tetranectin binds to plasminogen and to isolated kringle 4. May be involved in the packaging of molecules destined for exocytosis.,mass spectrometry: PubMed:10614823,online information:Tetranectin, similarity:Contains 1 C-type lectin domain.,subunit:Homotrimer.,tissue



specificity: Found in plasma.,


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200