



CD275 Polyclonal Antibody

Catalog No	YP-Ab-10741
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	ICOSLG B7H2 B7RP1 ICOSL KIAA0653
Immunogen	Synthetic peptide from human protein at AA range: 201-250
Specificity	The antibody detects endogenous CD275
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC-p 1:50-200, ELISA 1:10000-20000. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	B7H2;ICOSLG;B7 H2;B7RP-1;B7RP1;CD275;ICOSLG B7H2 B7RP1 ICOSL KIAA0653;ICOS ligand (B7 homolog 2) (B7-H2) (B7-like protein GI50) (B7-related protein 1) (B7RP-1) (CD antigen CD275);ICOS ligand (B7 homolog 2;B7-H2;B7-like protein GI50;B7-related protein 1;CD antigen CD275)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	Isoform 1 is widely expressed (brain, heart, kidney, liver, lung, pancreas, placenta, skeletal muscle, bone marrow, colon, ovary, prostate, testis, lymph nodes, leukocytes, spleen, thymus and tonsil), while isoform 2 is detected only in lymph nodes, leukocytes and spleen. Expressed on activated monocytes and dendritic cells.
Function	alternative products:Additional isoforms seem to exist,function:Ligand for the T-cell-specific cell surface receptor ICOS. Acts as a costimulatory signal for T-cell proliferation and cytokine secretion; induces also B-cell proliferation and differentiation into plasma cells. Could play an important role in mediating local tissue responses to inflammatory conditions, as well as in modulating the secondary immune response by co-stimulating memory T-cell function.,



induction: Constitutive expression is further enhanced by treatment with TNF- α in peripheral blood B-cells and monocytes, while it is decreased in dendritic cells., sequence caution: The sequence differs from that shown in position 300 onward for unknown reason., sequence caution: Translation N-terminally shortened., similarity: Belongs to the immunoglobulin superfamily. BTN/ MOG family., similarity: Contains 1 Ig-like C2-type (immunogl

Background

alternative products: Additional isoforms seem to exist, function: Ligand for the T-cell-specific cell surface receptor ICOS. Acts as a costimulatory signal for T-cell proliferation and cytokine secretion; induces also B-cell proliferation and differentiation into plasma cells. Could play an important role in mediating local tissue responses to inflammatory conditions, as well as in modulating the secondary immune response by co-stimulating memory T-cell function., induction: Constitutive expression is further enhanced by treatment with TNF- α in peripheral blood B-cells and monocytes, while it is decreased in dendritic cells., sequence caution: The sequence differs from that shown in position 300 onward for unknown reason., sequence caution: Translation N-terminally shortened., similarity: Belongs to the immunoglobulin superfamily. BTN/ MOG family., similarity: Contains 1 Ig-like C2-type (immunoglobulin-like) domain., similarity: Contains 1 Ig-like V-type (immunoglobulin-like) domain., tissue specificity: Isoform 1 is widely expressed (brain, heart, kidney, liver, lung, pancreas, placenta, skeletal muscle, bone marrow, colon, ovary, prostate, testis, lymph nodes, leukocytes, spleen, thymus and tonsil), while isoform 2 is detected only in lymph nodes, leukocytes and spleen. Expressed on activated monocytes and dendritic cells.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemical analysis of paraffin-embedded Human-tonsil, antibody was diluted at 1:100