



CTXL Polyclonal Antibody

Catalog No	YP-Ab-10755
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	VSIG2 CTH CTXL UNQ2770/PRO7154
Immunogen	Synthesized peptide derived from human CTXL. at AA range: 260-290
Specificity	This antibody detects endogenous levels of CTXL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000, ELISA 1:10000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	221413P1Rik;CT like protein;CTH;CTXL;VSIG2;VSIG2 CTH CTXL UNQ2770/PRO7154;V-set and immunoglobulin domain-containing protein 2 (Cortical thymocyte-like protein) (CT-like protein)
Observed Band	34kD
Calculated Molecular Weight	34kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Highly expressed in stomach, colon, prostate, trachea and thyroid glands and weakly in bladder and lung.
Function	similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain., similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,tissue specificity:Highly expressed in stomach, colon, prostate, trachea and thyroid glands and weakly in bladder and lung.,
Background	similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain., similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,tissue specificity:Highly expressed in stomach, colon, prostate, trachea and thyroid glands and weakly in bladder and lung.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of hela cells using primary antibody diluted at 1:1000(4°C overnight). Secondary antibody:Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25°C, 1 hour)