



FA2H Polyclonal Antibody

Catalog No	YP-Ab-10769
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	FA2H
Immunogen	Synthesized peptide derived from human FA2H. at AA range: 101-150
Specificity	FA2H Polyclonal Antibody detects endogenous levels of FA2H
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000, ELISA 1:10000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FAXDC1;Fatty acid hydroxylase domain-containing protein 1;Fatty acid alpha-hydroxylase;Fatty acid 2-hydroxylase;FAAH;FA2H;Fatty acid 2-hydroxylase (EC 1.-.-.-) (Fatty acid alpha-hydroxylase)
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Microsome membrane ; Multi-pass membrane protein .
Tissue Specificity	Detected in differentiating cultured keratinocytes (at protein level). Detected in epidermis and cultured keratinocytes (PubMed:17355976). Highly expressed in brain and colon. Detected at lower levels in testis, prostate, pancreas and kidney (PubMed:15337768).
Function	cofactor:Iron.,disease:Defects in FA2H are the cause of leukodystrophy dysmyelinating with spastic paraparesis with or without dystonia (DLDSP) [MIM:612443]. The disorder consists of a progressive neurologic disease manifested by spasticity, disordered tonicity of muscle, and white matter degeneration.,domain:The histidine box domains may contain the active site and/or be involved in metal ion binding.,function:Required for alpha-hydroxylation of free fatty acids and the formation of alpha-hydroxylated sphingolipids.,



induction:Up-regulated during keratinocyte differentiation.,similarity:Belongs to the SCS7 family.,similarity:Contains 1 cytochrome b5 heme-binding domain., tissue specificity:Detected in differentiating cultured keratinocytes (at protein level). Detected in epidermis and cultured keratinocytes. Highly expressed in brain and colon. Detected at lower levels in testis, prostate

Background

This gene encodes a protein that catalyzes the synthesis of 2-hydroxysphingolipids, a subset of sphingolipids that contain 2-hydroxy fatty acids. Sphingolipids play roles in many cellular processes and their structural diversity arises from modification of the hydrophobic ceramide moiety, such as by 2-hydroxylation of the N-acyl chain, and the existence of many different head groups. Mutations in this gene have been associated with leukodystrophy dysmyelinating with spastic paraparesis with or without dystonia.[provided by RefSeq, Mar 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of mouse-heart cells using primary antibody diluted at 1:2000(4°C overnight). Secondary antibody:Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25°C, 1 hour)