




Tmem72 Polyclonal Antibody

Catalog No	YP-Ab-10824
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	TMEM72 C10orf127 KSP37
Immunogen	Synthesized peptide derived from human Tmem72 Polyclonal
Specificity	This antibody detects endogenous levels of Tmem72.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000, ELISA 1:10000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	bA285G1.3;C1orf127;DKFZp686M585;KSP37;TMEM72;transmembrane protein 72;TMEM72 C10orf127 KSP37;Transmembrane protein 72 (Kidney-specific secretory protein of 37 kDa)
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	This gene encodes a transmembrane protein which may be expressed specifically in the kidney. [provided by RefSeq, Sep 2016],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



1 mouse-liver

2 CAC02

Western blot analysis of various lysate, antibody was diluted at 1000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000