



Serine/threonine-protein kinase US3 homolog

Polyclonal Antibody

Catalog No	YP-Ab-10874
Isotype	IgG
Reactivity	Human virus
Applications	ELISA
Gene Name	US3
Immunogen	Synthesized peptide derived from human Serine/threonine-protein kinase US3 homolog
Specificity	This antibody detects endogenous levels of human Serine/threonine-protein kinase US3 homolog
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	ELISA(peptide) 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	US3;Serine/threonine-protein kinase US3 homolog
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	P09251 tegument protein; phosphorylates nuclear egress lamina protein; mediates phosphorylation of HDAC1 and HDAC2 and other cellular and viral proteins(ORF66) Human herpesvirus 3 catalytic activity:ATP + a protein = ADP + a phosphoprotein.,function:Phosphorylates IE62.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family.,similarity:Contains 1 protein kinase domain.,

Background

P09251 tegument protein; phosphorylates nuclear egress lamina protein; mediates phosphorylation of HDAC1 and HDAC2 and other cellular and viral proteins(ORF66) Human herpesvirus 3 catalytic activity:ATP + a protein = ADP + a phosphoprotein.,function:Phosphorylates IE62.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family.,similarity:Contains 1 protein kinase domain.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images