



Serine/threonine-protein kinase UL13 homolog

Polyclonal Antibody

Catalog No	YP-Ab-10875
Isotype	IgG
Reactivity	Human virus
Applications	ELISA
Gene Name	UL13
Immunogen	Synthesized peptide derived from human Serine/threonine-protein kinase UL13 homolog
Specificity	This antibody detects endogenous levels of human Serine/threonine-protein kinase UL13 homolog
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	ELISA(peptide) 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UL13;Serine/threonine-protein kinase UL13 homolog
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images