



Probable ganciclovir kinase U69 Polyclonal Antibody

Catalog No	YP-Ab-10882
Isotype	IgG
Reactivity	Human virus
Applications	ELISA
Gene Name	U69
Immunogen	Synthesized peptide derived from human Probable ganciclovir kinase U69
Specificity	This antibody detects endogenous levels of human Probable ganciclovir kinase U69
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	ELISA(peptide) 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Probable ganciclovir kinase U69;U69
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	P52344 tegument serine/threonine protein kinase(U69) Human herpesvirus 7 function:Phosphorylates the antiviral nucleoside analog ganciclovir., similarity:Belongs to the protein kinase superfamily. Tyr protein kinase family. HCMV ganciclovir subfamily.,
Background	P52344 tegument serine/threonine protein kinase(U69) Human herpesvirus 7 function:Phosphorylates the antiviral nucleoside analog ganciclovir., similarity:Belongs to the protein kinase superfamily. Tyr protein kinase family. HCMV ganciclovir subfamily.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images