




ELOF1 Polyclonal Antibody

Catalog No	YP-Ab-10908
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	IHC; ELISA
Gene Name	ELOF1
Immunogen	Synthesized peptide derived from human ELOF1 AA range: 30-110
Specificity	This antibody detects endogenous levels of human ELOF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ELF1; ELOF1; Transcription elongation factor 1 homolog
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Muscle,
Function	function: Transcription elongation factor implicated in the maintenance of proper chromatin structure in actively transcribed regions., similarity: Belongs to the ELOF1 family.,
Background	function: Transcription elongation factor implicated in the maintenance of proper chromatin structure in actively transcribed regions., similarity: Belongs to the ELOF1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).