




TM4SF19 Polyclonal Antibody

Catalog No	YP-Ab-10928
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB; ELISA
Gene Name	TM4SF19 OCTM4
Immunogen	Synthesized peptide derived from human TM4SF19. AA range:140-160
Specificity	This antibody detects endogenous levels of human TM4SF19
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000, ELISA(peptide)1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	OCTM4;TM4SF19;TM4SF19 OCTM4;Transmembrane 4 L6 family member 19 (Osteoclast maturation-associated gene 4 protein;Tetraspan membrane protein OCTM4)
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Skin,
Function	similarity:Belongs to the L6 tetraspanin family.,
Background	similarity:Belongs to the L6 tetraspanin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of various lysates, primary antibody was diluted at 1:1000, 4° over night, secondary antibody(cat: RS23920)was diluted at 1:10000, 37° 1hour.