




IER5 rabbit pAb

Catalog No	YP-Ab-10933
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	IER5 PP4583 SBB148
Immunogen	Synthesized peptide derived from human IER5 AA range: 28-78
Specificity	This antibody detects endogenous levels of IER5 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SBB148;IER5;IER5 PP4583 SBB148
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Cytoplasm . Predominantly cytoplasmic (PubMed:25816751). Translocated in the nucleus during heat shock (PubMed:25816751). .
Tissue Specificity	Expressed in acute myeloid leukemia (AML) cells.
Function	similarity:Belongs to the IER family.,
Background	This gene encodes a protein that is similar to other immediate early response proteins. In the mouse, a similar gene may play an important role in mediating the cellular response to mitogenic signals. Studies in rats found the expression of a similar gene to be increased after waking and sleep deprivation. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from AD293 cells, primary antibody was diluted at 1:1000, 4° over night