




F198B rabbit pAb

Catalog No	YP-Ab-10934
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	FAM198B C4orf18 ENED AD021 UNQ2512/PRO6001
Immunogen	Synthesized peptide derived from human F198B AA range: 372-422
Specificity	This antibody detects endogenous levels of F198B at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F198B;FAM198B C4orf18 ENED AD021 UNQ2512/PRO6001
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Single-pass type II membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4°over night