




CA159 rabbit pAb

Catalog No	YP-Ab-10961
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C1orf159 UNQ2998/PRO9739
Immunogen	Synthesized peptide derived from human CA159 AA range: 245-295
Specificity	This antibody detects endogenous levels of CA159 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CA159;C1orf159 UNQ2998/PRO9739
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4°over night