




F207A rabbit pAb

Catalog No	YP-Ab-10966
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM207A C21orf70 PRED56
Immunogen	Synthesized peptide derived from human F207A AA range: 29-79
Specificity	This antibody detects endogenous levels of F207A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F207A;FAM207A C21orf70 PRED56
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Not detected in any tested tissue.
Function	tissue specificity:Not detected in any tested tissue.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4°over night