




KCTD6 rabbit pAb


Catalog No	YP-Ab-10967
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;IHC
Gene Name	KCTD6
Immunogen	Synthesized peptide derived from human KCTD6 AA range: 118-168
Specificity	This antibody detects endogenous levels of KCTD6 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KCTD6
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, myofibril, sarcomere, M line . Colocalizes with ANK1 isoform Mu17 at the M line in differentiated skeletal muscle cells and heart. .
Tissue Specificity	Highly expressed in cerebellum and brain. Expression is down-regulated in medulloblastoma.
Function	similarity:Contains 1 BTB (POZ) domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).