




LNP1 rabbit pAb

Catalog No	YP-Ab-10974
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LNP1 NP3
Immunogen	Synthesized peptide derived from human LNP1 AA range: 11-61
Specificity	This antibody detects endogenous levels of LNP1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LNP1;LNP1 NP3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	disease:A chromosomal aberration involving LNP1 is found in a form of T-cell acute lymphoblastic leukemia (T-ALL). Translocation t(3;11)(q12.2;p15.4) with NUP98.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4°over night