




GLYR1 rabbit pAb

Catalog No	YP-Ab-10980
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	GLYR1 HIBDL NP60
Immunogen	Synthesized peptide derived from human GLYR1 AA range: 130-180
Specificity	This antibody detects endogenous levels of GLYR1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	3-hydroxyisobutyrate dehydrogenase-like protein; Cytokine-like nuclear factor N-PAC; GLYR1; GLYR1 HIBDL NP60
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Chromosome .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from CACO2 cells,
primary antibody was diluted at 1:1000, 4°over night