




CG060 rabbit pAb

Catalog No	YP-Ab-11016
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C7orf60
Immunogen	Synthesized peptide derived from human CG060 AA range: 108-158
Specificity	This antibody detects endogenous levels of CG060 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CG060;C7orf60
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleolus,
Tissue Specificity	
Function	similarity:Belongs to the UPF0532 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4°over night