




LUR1L rabbit pAb

Catalog No	YP-Ab-11020
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	LURAP1L C9orf150 HYST0841
Immunogen	Synthesized peptide derived from human LUR1L AA range: 142-192
Specificity	This antibody detects endogenous levels of LUR1L at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LUR1L; LURAP1L C9orf150 HYST0841
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	polymorphism: The poly-Gly region is polymorphic and the number of Gly varies in the population (from 9 to 12).,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night