




ANKZ1 rabbit pAb

	YP-Ab-11021
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	ANKZF1 ZNF744
Immunogen	Synthesized peptide derived from human ANKZ1 AA range: 478-528
Specificity	This antibody detects endogenous levels of ANKZ1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANKZ1; ANKZF1 ZNF744
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Translocates to the mitochondria upon exposure to hydrogen peroxide. .
Tissue Specificity	
Function	similarity: Contains 1 C2H2-type zinc finger., similarity: Contains 2 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night