




CG062 rabbit pAb

Catalog No	YP-Ab-11023
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	C7orf62
Immunogen	Synthesized peptide derived from human CG062 AA range: 59-109
Specificity	This antibody detects endogenous levels of CG062 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CG062; C7orf62
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Testis-specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night