




GUC2A rabbit pAb

Catalog No	YP-Ab-11050
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	GUCA2A GUCA2
Immunogen	Synthesized peptide derived from human GUC2A AA range: 52-102
Specificity	This antibody detects endogenous levels of GUC2A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANP A; ANPa; ANPR A; ANPRA; GC A; Guanylate cyclase A; GUC2A; GUCY2A; NPR1; GUCA2A GUCA2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted.
Tissue Specificity	Highly expressed in ileum and colon. Found in plasma.
Function	function: Endogenous activator of intestinal guanylate cyclase. It stimulates this enzyme through the same receptor binding region as the heat-stable enterotoxins., similarity: Belongs to the guanylin family., tissue specificity: Highly expressed in ileum and colon. Found in plasma.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HeLa cells,
primary antibody was diluted at 1:1000, 4° over night