




NUDC2 rabbit pAb

Catalog No	YP-Ab-11053
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	NUDCD2
Immunogen	Synthesized peptide derived from human NUDC2 AA range: 90-140
Specificity	This antibody detects endogenous levels of NUDC2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NUDC2;NUDCD2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Chromosome, centromere, kinetochore . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle pole . Associates with centrosomes in interphase and to spindle poles and kinetochores during mitosis.
Tissue Specificity	
Function	similarity:Contains 1 CS domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4°over night