




PPR3E rabbit pAb

Catalog No	YP-Ab-11074
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	PPP1R3E
Immunogen	Synthesized peptide derived from human PPR3E AA range: 210-260
Specificity	This antibody detects endogenous levels of PPR3E at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PPR3E; PPP1R3E
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in skeletal muscle and heart with barely detectable levels in liver.
Function	domain: The CBM21 domain is known to be involved in the localization to glycogen and is characteristic of some regulatory subunit of phosphatase complexes., function: Acts as a glycogen-targeting subunit for PP1. PP1 is involved in glycogen metabolism and contributes to the activation of glycogen synthase leading to an increase in glycogen synthesis., similarity: Contains 1 CBM21 (carbohydrate binding type-21) domain., tissue specificity: Expressed in skeletal muscle and heart with barely detectable levels in liver.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night