




FA54A rabbit pAb

Catalog No	YP-Ab-11086
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	FAM54A DUFD1
Immunogen	Synthesized peptide derived from human FA54A AA range: 303-353
Specificity	This antibody detects endogenous levels of FA54A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA54A; FAM54A DUFD1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion . Associated with membranes. .
Tissue Specificity	
Function	similarity: Belongs to the MTFR1/FAM54 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night