




SELV rabbit pAb

Catalog No	YP-Ab-11087
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SELV
Immunogen	Synthesized peptide derived from human SELV AA range: 250-300
Specificity	This antibody detects endogenous levels of SELV at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SELV
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Testis specific.
Function	function:May be involved in a redox-related process .,similarity:Belongs to the SelWTH family.,tissue specificity:Testis specific.,
Background	This gene encodes a selenoprotein, which contains a selenocysteine (Sec) residue at its active site. The selenocysteine is encoded by the UGA codon that normally signals translation termination. The 3' UTR of selenoprotein genes have a common stem-loop structure, the sec insertion sequence (SECIS), that is necessary for the recognition of UGA as a Sec codon rather than as a stop signal. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night