




PAQR9 rabbit pAb

Catalog No	YP-Ab-11105
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PAQR9
Immunogen	Synthesized peptide derived from human PAQR9 AA range: 68-118
Specificity	This antibody detects endogenous levels of PAQR9 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PAQR9
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Expression levels vary widely in a range of tissues (PubMed:16044242). Expressed in the brain, at high level in the pituitary gland and also in hypothalamus, limbic system, caudate nucleus accumbens, pons and olfactory bulb (PubMed:23161870).
Function	similarity:Belongs to the ADIPOR family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night