




IGHG4 rabbit pAb

Catalog No	YP-Ab-11106
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	IGHG4
Immunogen	Synthesized peptide derived from human IGHG4 AA range: 178-228
Specificity	This antibody detects endogenous levels of IGHG4 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IGHG4;Ig gamma-4 chain C region;Human IgG4/Immunoglobulin heavy constant gamma 4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted . Cell membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night