




NPPL3 rabbit pAb

Catalog No	YP-Ab-11111
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NPIPL3 KIAA0220
Immunogen	Synthesized peptide derived from human NPPL3 AA range: 831-881
Specificity	This antibody detects endogenous levels of NPPL3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NPPL3;NPIPL3 KIAA0220
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	caution:Could be the product of a pseudogene.,similarity:Belongs to the NPIP family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4°over night