




FDXA1 rabbit pAb

Catalog No	YP-Ab-11143
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FDXACB1
Immunogen	Synthesized peptide derived from human FDXA1 AA range: 358-408
Specificity	This antibody detects endogenous levels of FDXA1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FDXA1;FDXACB1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoplasm,
Tissue Specificity	
Function	similarity:Contains 1 FDX-ACB domain.,
Background	This gene encodes a protein which contains a ferredoxin-fold anticodon-binding domain which is contained in a subset of phenylalanyl tRNA synthetases. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC-12 cells,
primary antibody was diluted at 1:1000, 4° over night