




ACPL2 rabbit pAb

Catalog No	YP-Ab-11144
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	ACPL2 UNQ370/PRO706
Immunogen	Synthesized peptide derived from human ACPL2 AA range: 201-251
Specificity	This antibody detects endogenous levels of ACPL2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACPL2; ACPL2 UNQ370/PRO706
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Single-pass type II membrane protein . Colocalizes to Golgi apparatus in a B3GAT3-dependent manner. .
Tissue Specificity	Widely expressed. Strongly expressed in spleen, fetal liver, moderately in placenta, pancreas, kidney, thymus and colon.
Function	catalytic activity: A phosphate monoester + H ₂ O = an alcohol + phosphate., similarity: Belongs to the histidine acid phosphatase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night