




MZT1 rabbit pAb

Catalog No	YP-Ab-11145
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	MZT1 C13orf37 MOZART1
Immunogen	Synthesized peptide derived from human MZT1 AA range: 19-69
Specificity	This antibody detects endogenous levels of MZT1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C13orf37;MOZART1;MZT1;MZT1 C13orf37 MOZART1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle .
Tissue Specificity	
Function	similarity:Belongs to the UPF0582 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).