



F165B rabbit pAb

Catalog No	YP-Ab-11153
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	FAM165B C21orf51
Immunogen	Synthesized peptide derived from human F165B AA range: 7-57
Specificity	This antibody detects endogenous levels of F165B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F165B;FAM165B C21orf51
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	Expressed in heart, spleen, liver, stomach, muscle, lung, testis, skin, PBL and bone marrow.
Function	similarity:Belongs to the FAM165 family.,tissue specificity:Expressed in heart, spleen, liver, stomach, muscle, lung, testis, skin, PBL and bone marrow.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images