




SPAC7 rabbit pAb


Catalog No	YP-Ab-11175
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;ELISA;IHC
Gene Name	SPACA7 C13orf28
Immunogen	Synthesized peptide derived from human SPAC7 AA range: 115-165
Specificity	This antibody detects endogenous levels of SPAC7 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPAC7;SPACA7 C13orf28
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted . Cytoplasmic vesicle, secretory vesicle, acrosome . Cytoplasmic vesicle, secretory vesicle, acrosome lumen . Localized in perinuclear pro-acrosomal granules in round spermatides. Localized between the inner and outer acrosomal membranes (matrix or lumen) in spermatozoa. Secreted during acrosome exocytosis. .
Tissue Specificity	Expressed in spermatozoa (at protein level) (PubMed:22495889).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from VEC cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human oophoroma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).