




FA73A rabbit pAb

Catalog No	YP-Ab-11177
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM73A
Immunogen	Synthesized peptide derived from human FA73A AA range: 36-86
Specificity	This antibody detects endogenous levels of FA73A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA73A;FAM73A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion outer membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	caution:It is uncertain whether Met-1 or Met-33 is the initiator.,similarity:Belongs to the FAM73 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night