




T53G1 rabbit pAb

Catalog No	YP-Ab-11185
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF
Gene Name	TP53TG1 P53TG1 TP53AP1
Immunogen	Synthesized peptide derived from human T53G1 AA range: 19-69
Specificity	This antibody detects endogenous levels of T53G1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T53G1;TP53TG1 P53TG1 TP53AP1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	induction:Up-regulated by ultraviolet irradiation and exposure to bleomycin or cisplatin in a TP53/p53-dependent manner.,tissue specificity:Isoform 1, isoform 2, isoform 3 and isoform 4 are ubiquitously expressed.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).