



HEM6 rabbit pAb

Catalog No	YP-Ab-11189
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	CPOX CPO CPX
Immunogen	Synthesized peptide derived from human HEM6 AA range: 265-315
Specificity	This antibody detects endogenous levels of HEM6 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HEM6; CPOX CPO CPX
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion intermembrane space.
Tissue Specificity	
Function	catalytic activity: Coproporphyrinogen-III + O(2) + 2 H(+) = protoporphyrinogen-IX + 2 CO(2) + 2 H(2)O., disease: Defects in CPOX are the cause of hereditary coproporphyria (HCP) [MIM:121300]. HCP is an acute hepatic porphyria and an autosomal dominant disease characterized by neuropsychiatric disturbances and skin photosensitivity. Biochemically, there is an overexcretion of coproporphyrin III in the urine and in the feces. HCP is clinically characterized by attacks of abdominal pain, neurological disturbances, and psychiatric symptoms. The symptoms are generally manifested with rapid onset, and can be precipitated by drugs, alcohol, caloric deprivation, infection, endocrine factors or stress. A severe variant form is harderoporphyria, which is characterized by earlier onset attacks, massive excretion of harderoporphyria in the feces, and a marked decrease of coproporphyrinogen IX oxidase



Background

The protein encoded by this gene is the sixth enzyme of the heme biosynthetic pathway. The encoded enzyme is soluble and found in the intermembrane space of mitochondria. This enzyme catalyzes the stepwise oxidative decarboxylation of coproporphyrinogen III to protoporphyrinogen IX, a precursor of heme. Defects in this gene are a cause of hereditary coproporphyria (HCP) .[provided by RefSeq, Oct 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from PC-12 cells, primary antibody was diluted at 1:1000, 4° over night