




WDR54 rabbit pAb

Catalog No	YP-Ab-11191
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	WDR54
Immunogen	Synthesized peptide derived from human WDR54 AA range: 99-149
Specificity	This antibody detects endogenous levels of WDR54 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	WDR45;WDRX1;WD repeat-containing protein 45;WD repeat domain phosphoinositide-interacting protein 4;WD repeat;WDR54
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Vesicle . Aggregates in vesicles when cross-linked. .
Tissue Specificity	Expressed in epithelial cells (at protein level) (PubMed:29987896). Isoform 3 expression is highly increased in colorectal cancer cells (PubMed:29987896).
Function	similarity:Contains 4 WD repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4°over night