




BN3D2 rabbit pAb

Catalog No	YP-Ab-11193
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	BCDIN3D
Immunogen	Synthesized peptide derived from human BN3D2 AA range: 192-242
Specificity	This antibody detects endogenous levels of BN3D2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BN3D2;BCDIN3D
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	function:Probable methyltransferase.,similarity:Belongs to the methyltransferase superfamily.,similarity:Contains 1 BIN3 domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Jarkat cells,
primary antibody was diluted at 1:1000, 4°over night