




CU024 rabbit pAb

Catalog No	YP-Ab-11203
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LINC00114 C21orf24 NCRNA00114
Immunogen	Synthesized peptide derived from human CU024 AA range: 4-54
Specificity	This antibody detects endogenous levels of CU024 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CU024;LINC00114 C21orf24 NCRNA00114
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in a range of cell lines, including B-cell lymphoma and prostate.
Function	caution:Product of a dubious CDS prediction. Probable non-coding RNA. At least 6 alternative mRNA transcripts exist, with consensus splice donor and acceptor sites, but no coding potential nor murine ortholog. This protein sequence corresponds to the so-called isoforms C21orf24-6 and C21orf24-7., tissue specificity:Expressed in a range of cell lines, including B-cell lymphoma and prostate.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night