



ZBED6 rabbit pAb

Catalog No	YP-Ab-11208
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZBED6
Immunogen	Synthesized peptide derived from human ZBED6 AA range :130-190
Specificity	This antibody detects endogenous levels of ZBED6 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZBED6
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Nucleus, nucleolus . Cytoplasm . Located predominantly in the nucleolus but is also dispersed throughout the nucleus (By similarity). Mainly cytoplasmic in insulin- and glucagon-positive islet pancreatic cells, however nuclear localization is evidence in some pancreatic endocrine cells (PubMed:24043816). .
Tissue Specificity	Expressed in pancreatic islet cells (at protein level).
Function	
Background	The protein encoded by this transposon-derived intronless gene is a transcriptional repressor that binds to the consensus sequence 5'-GCTCGC-3'. The encoded protein has been shown to repress IGF2 transcription. This gene is located within the first intron of the ZC3H11A gene. [provided by RefSeq, Jul 2016],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HEK293 cells, primary antibody was diluted at 1:1000, 4° over night