




PDZ1P rabbit pAb

Catalog No	YP-Ab-11221
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	PDZK1P1; PDZK1P2
Immunogen	Synthesized peptide derived from human PDZ1P AA range: 297-347
Specificity	This antibody detects endogenous levels of PDZ1P at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PDZ1P;PDZK1P1;PDZK1P2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4°over night