




NO66 rabbit pAb

Catalog No	YP-Ab-11227
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	NO66 C14orf169 MAPJD
Immunogen	Synthesized peptide derived from human NO66 AA range: 533-583
Specificity	This antibody detects endogenous levels of NO66 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C14orf169; hsNO66; MAPJD; NO66; NO66/C14orf169; Nucleolar protein 66; NO66 C14orf169 MAPJD
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus. Nucleus, nucleoplasm. Granular part of nucleoli. Nucleoplasm, nucleoplasmic foci, some of them associated with nucleoli.
Tissue Specificity	Widely expressed. Overexpressed in lung carcinomas.
Function	function: May play a role in ribosome biogenesis and in the replication or remodeling of certain heterochromatic region., similarity: Belongs to the MINA53/NO66 family., similarity: Contains 1 JmjC domain., subcellular location: Granular part of nucleoli. Nucleoplasm, nucleoplasmic foci, some of them associated with nucleoli., tissue specificity: Ubiquitous.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night